

For Research Use Only

Mouse Il1r2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85071-3-PBS



Basic Information

Catalog Number:

85071-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_010555.4

GeneID (NCBI):

16178

UNIPROT ID:

P27931

Full Name:

interleukin 1 receptor, type II

Calculated MW:

46 kDa

Purification Method:

Protein A purification

CloneNo.:

242602F1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01814-2, MOUSE IL-1R2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85071-1-PBS. Detection antibody: 85071-3-PBS. Standard: Eg2956. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01814-1, MOUSE IL-1R2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85071-2-PBS. Detection antibody: 85071-3-PBS. Standard: Eg2956. Range: 0.781-100 ng/mL