

For Research Use Only

Mouse Il1r2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85071-5-PBS



Basic Information

Catalog Number:	85071-5-PBS	GenBank Accession Number:	NM_010555.4	Purification Method:	Protein A purification
Concentration:	1 mg/ml	GeneID (NCBI):	16178	CloneNo.:	242602H6
Source:	Rabbit	UNIPROT ID:	P27931		
Isotype:	IgG	Full Name:	interleukin 1 receptor, type II		
Immunogen Catalog Number:	EG2956	Calculated MW:	46 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean IL-1R2 concentration was determined to be 208.4 ng/mL in mouse serum.



Cytometric bead array standard curve of MP01814-4, MOUSE IL-1R2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85071-2-PBS. Detection antibody: 85071-5-PBS. Standard: Eg2956. Range: 0.781-100 ng/mL.