

For Research Use Only

Mouse Slamf1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85073-4-PBS



Basic Information

Catalog Number:

85073-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_013730.4

GeneID (NCBI):

27218

UNIPROT ID:

Q9QUM4-1

Full Name:

signaling lymphocytic activation
molecule family member 1

Calculated MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

242631G7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01792-2, MOUSE SLAM/CD150 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85073-3-PBS. Detection antibody: 85073-4-PBS. Standard: Eg2787. Range: 0.625-40 ng/mL