

For Research Use Only

# FLT4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85076-1-PBS



## Basic Information

Catalog Number:

85076-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_182925

GeneID (NCBI):

2324

UNIPROT ID:

P35916-1

Full Name:

fms-related tyrosine kinase 4

Calculated MW:

153kd

Purification Method:

Protein A purification

CloneNo.:

242430C2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00009-4, FLT4/VEGFR3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85076-1-PBS. Detection antibody: 85076-3-PBS. Standard: Eg2494. Range: 0.781-100 ng/mL.