

For Research Use Only

GORASP2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85077-1-PBS



Basic Information

Catalog Number:

85077-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0950

GenBank Accession Number:

BC007770

GeneID (NCBI):

26003

UNIPROT ID:

Q9H8Y8

Full Name:

golgi reassembly stacking protein 2,

55kDa

Calculated MW:

55 kDa

Purification Method:

Protein A purification

CloneNo.:

242700B5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01795-2, GORASP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85077-1-PBS. Detection antibody: 85077-2-PBS. Standard: Ag0950. Range: 0.781-100 ng/mL.



Biolayer interferometry (BLI) kinetic assays of 85077-1-RR against Human GORASP2 were performed. The affinity constant is 0.246 nM.