

For Research Use Only

ZP3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85084-2-PBS



Basic Information

Catalog Number:

85084-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15791

GenBank Accession Number:

BC113949

GeneID (NCBI):

7784

UNIPROT ID:

P21754

Full Name:

zona pellucida glycoprotein 3 (sperm receptor)

Calculated MW:

373 aa, 41 kDa

Purification Method:

Protein A purification

CloneNo.:

242729B9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01801-2, ZP3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85084-1-PBS. Detection antibody: 85084-2-PBS. Standard: Ag15791. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01801-1, ZP3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85084-3-PBS. Detection antibody: 85084-2-PBS. Standard: Ag15791. Range: 0.781-100 ng/mL.