

For Research Use Only

ASIC1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85089-2-PBS



Basic Information

Catalog Number:

85089-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25941

GenBank Accession Number:

NM_001095

GeneID (NCBI):

41

UNIPROT ID:

P78348

Full Name:

amiloride-sensitive cation channel 2,
neuronal

Calculated MW:

60 kDa

Purification Method:

Protein A purification

CloneNo.:

242749F9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01789-2, ASIC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85089-1-PBS. Detection antibody: 85089-2-PBS. Standard: Ag25941. Range: 0.313-40 ng/mL