

For Research Use Only

Galectin 1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 85091-2-PBS



Basic Information

Catalog Number:

85091-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG7626

GenBank Accession Number:

BC020675

GeneID (NCBI):

3956

UNIPROT ID:

P09382

Full Name:

lectin, galactoside-binding, soluble, 1

Calculated MW:

15 kDa

Purification Method:

Protein A purification

CloneNo.:

260418C6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Galectin 1 concentration was determined to be 22.6 ng/mL in Human plasma 1, 19.5 ng/mL in Human plasma 2.



Cytometric bead array standard curve of MP03369-1, Galectin 1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85091-3-PBS. Detection antibody: 85091-2-PBS. Standard: Eg7626. Range: 1.562-200 ng/mL.