

For Research Use Only

THRAP3 Recombinant antibody, PBS Only

Catalog Number: 85093-3-PBS



Basic Information

Catalog Number:

85093-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_005119

GeneID (NCBI):

9967

UNIPROT ID:

Q9Y2W1

Full Name:

thyroid hormone receptor associated protein 3

Calculated MW:

109 kDa

Observed MW:

109-130 kDa

Purification Method:

Protein A purification

CloneNo.:

242268A9

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85093-3-RR (THRAP3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85093-3-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 85093-3-RR (THRAP3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85093-3-PBS in a different storage buffer formulation.