

For Research Use Only

ND4 Recombinant antibody

Catalog Number: 85098-1-RR



Basic Information

Catalog Number:

85098-1-RR

Concentration:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

YP_003024035

GeneID (NCBI):

4538

UNIPROT ID:

P03905

Full Name:

NADH dehydrogenase, subunit 4
(complex I)

Calculated MW:

52 kDa

Observed MW:

43-52 kDa

Purification Method:

Protein A purification

CloneNo.:

242422E2

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells, HepG2 cells, SGC-7901 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85098-1-RR (ND4 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.