

For Research Use Only

SUN2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85103-1-PBS



Basic Information

Catalog Number: 85103-1-PBS	GenBank Accession Number: BC111717	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 25777	CloneNo.: 242533F8
Source: Rabbit	UNIPROT ID: Q9UH99	
Isotype: IgG	Full Name: unc-84 homolog B (C. elegans)	
Immunogen Catalog Number: AG16077	Calculated MW: 717 aa, 80 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01849-1, SUN2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85103-3-PBS. Detection antibody: 85103-1-PBS. Standard: Ag16077. Range: 0.625-80 ng/mL