

For Research Use Only

Caspase 10 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85104-1-PBS



Basic Information

Catalog Number:

85104-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC042844

GeneID (NCBI):

843

UNIPROT ID:

Q92851

Full Name:

caspase 10, apoptosis-related
cysteine peptidase

Calculated MW:

59 kDa

Purification Method:

Protein A purification

CloneNo.:

242468B5

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01855-1, Human Caspase 10 Recombinant Matched Antibody Pair - PBS only. 85104-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard SY02408. 85104-1-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.



The mean Caspase 10 concentration was determined to be 4.6 ng/mL in HeLa cell extract based on a 1.0 mg/mL extract load.