

For Research Use Only

KIRREL3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85106-1-PBS



Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
85106-1-PBS	BC101801	Protein A purification
Concentration:	GeneID (NCBI):	CloneNo.:
1000 ug/ml	84623	242816F11
Source:	UNIPROT ID:	
Rabbit	Q8IZU9	
Isotype:	Full Name:	
IgG	kin of IRRE like 3 (Drosophila)	
Immunogen Catalog Number:	Calculated MW:	
AG12573	778 aa, 85 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01859-1, KIRREL3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85106-2-PBS. Detection antibody: 85106-1-PBS. Standard: Ag12573. Range: 0.781-100 ng/mL