

For Research Use Only

SYTL5 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85108-1-PBS



Basic Information

Catalog Number:

85108-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21698

GenBank Accession Number:

BC131585

GeneID (NCBI):

94122

UNIPROT ID:

Q8TDW5

Full Name:

synaptotagmin-like 5

Calculated MW:

730 aa, 81 kDa

Purification Method:

Protein A purification

CloneNo.:

242838D3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01854-2, SYTL5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85108-1-PBS. Detection antibody: 85108-4-PBS. Standard: Ag21698. Range: 1.563-200 ng/mL.