

For Research Use Only

MRPL11 Recombinant antibody, PBS Only

Catalog Number: 85109-5-PBS



Basic Information

Catalog Number:

85109-5-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7899

GenBank Accession Number:

BC005002

GeneID (NCBI):

65003

UNIPROT ID:

Q9Y3B7

Full Name:

mitochondrial ribosomal protein L11

Calculated MW:

21 kDa

Observed MW:

21 kDa

Purification Method:

Protein A purification

CloneNo.:

242768G7

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85109-5-RR (MRPL11 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85109-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 85109-5-RR against Human MRPL11 were performed. The affinity constant is 0.324 nM.