

For Research Use Only

SEC13 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85110-1-PBS



Basic Information

Catalog Number:

85110-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7610

GenBank Accession Number:

BC002634

GeneID (NCBI):

6396

UNIPROT ID:

P55735

Full Name:

SEC13 homolog (S. cerevisiae)

Calculated MW:

36 kDa

Purification Method:

Protein A purification

CloneNo.:

242775E3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01867-2, SEC13 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85110-1-PBS. Detection antibody: 85110-3-PBS. Standard: Ag7610. Range: 0.781-100 ng/mL.