

For Research Use Only

PR Recombinant antibody, PBS Only (Detector)

Catalog Number: 85112-1-PBS



Basic Information

Catalog Number:

85112-1-PBS

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22986

GenBank Accession Number:

NM_000926

GeneID (NCBI):

5241

UNIPROT ID:

P06401

Full Name:

progesterone receptor

Calculated MW:

933 aa, 99 kDa

Purification Method:

Protein A purification

CloneNo.:

242686B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01857-1, Progesterone receptor Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85112-3-PBS. Detection antibody: 85112-1-PBS. Standard: Ag22986. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01857-2, Progesterone receptor Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85112-2-PBS. Detection antibody: 85112-1-PBS. Standard: Ag22986. Range: 1.25-80 ng/mL