

For Research Use Only

HCLS1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85114-3-PBS



Basic Information

Catalog Number:

85114-3-PBS

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8645

GenBank Accession Number:

BC016758

GeneID (NCBI):

3059

UNIPROT ID:

P14317

Full Name:

hematopoietic cell-specific Lyn
substrate 1

Calculated MW:

486 aa, 54 kDa

Purification Method:

Protein A purification

CloneNo.:

242825G11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01860-2, HCLS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85114-1-PBS. Detection antibody: 85114-3-PBS. Standard: Ag8645. Range: 0.781-100 ng/mL.