

For Research Use Only

Galectin 3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85115-1-PBS



Basic Information

Catalog Number:

85115-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0227

GenBank Accession Number:

BC001120

GeneID (NCBI):

3958

UNIPROT ID:

P17931

Full Name:

lectin, galactoside-binding, soluble, 3

Calculated MW:

26 kDa

Purification Method:

Protein A purification

CloneNo.:

242291A8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01833-1, Galectin-3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85115-1-PBS. Detection antibody: 85115-2-PBS. Standard: Eg0227. Range: 0.625-40 ng/mL