

For Research Use Only

# TPX2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85118-2-PBS



## Basic Information

Catalog Number:

85118-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2334

GenBank Accession Number:

BC020207

GeneID (NCBI):

22974

UNIPROT ID:

Q9ULW0

Full Name:

TPX2, microtubule-associated,  
homolog (Xenopus laevis)

Calculated MW:

747 aa, 86 kDa

Purification Method:

Protein A purification

CloneNo.:

242528C11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01821-1, TPX2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85118-2-PBS. Detection antibody: 85118-3-PBS. Standard: Ag2334. Range: 0.625-40 ng/mL