

For Research Use Only

FOXA1 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 85125-3-PBS



Basic Information

Catalog Number: 85125-3-PBS	GenBank Accession Number: BC033890	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 3169	CloneNo.: 242658C6
Source: Rabbit	UNIPROT ID: P55317	
Isotype: IgG	Full Name: forkhead box A1	
Immunogen Catalog Number: AG14243	Calculated MW: 473 aa, 49 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01842-2, FOXA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85125-1-PBS. Detection antibody: 85125-3-PBS. Standard: Ag14243. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP01842-1, FOXA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85125-3-PBS. Detection antibody: 85125-2-PBS. Standard: Ag14243. Range: 0.156-20 ng/mL