

For Research Use Only

Mouse CD36/SR-B3 Recombinant monoclonal antibody, PBS Only

Catalog Number: 85127-4-PBS



Basic Information

Catalog Number:

85127-4-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3464

GenBank Accession Number:

NM_007643.4

GeneID (NCBI):

12491

UNIPROT ID:

Q08857

Full Name:

CD36 antigen

Calculated MW:

53 kDa

Purification Method:

Protein A purification

CloneNo.:

242662F7

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

