

For Research Use Only

# Foxp3 Recombinant antibody

Catalog Number: 85128-4-RR



## Basic Information

Catalog Number:

85128-4-RR

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

20371

Full Name:

forkhead box P3

Observed MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

242674B10

Recommended Dilutions:

WB 1:5000-1:50000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : MJ cells, 4T1 cells, RAW 264.7 cells.

## Background Information

FOXP3, also named as IPEX, JM2 and Scurfin, is a transcription factor. FOXP3 plays a critical role in the control of immune response. Defects in FOXP3 are the cause of immunodeficiency polyendocrinopathy, enteropathy, X-linked syndrome (IPEX) which also known as X-Linked autoimmunity-immunodeficiency syndrome.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85128-4-RR (Foxp3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.