

For Research Use Only

ASCC3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85130-4-PBS



Basic Information

Catalog Number:

85130-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11848

GenBank Accession Number:

BC050681

GeneID (NCBI):

10973

UNIPROT ID:

Q8N3C0

Full Name:

activating signal cointegrator 1
complex subunit 3

Calculated MW:

111 aa, 13 kDa, 251 kDa

Purification Method:

Protein A purification

CloneNo.:

242755F1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01828-2, ASCC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85130-4-PBS. Detection antibody: 85130-3-PBS. Standard: Ag11848. Range: 0.781-100 ng/mL.