

For Research Use Only

# ASCC3 Recombinant antibody, PBS Only

Catalog Number: 85130-5-PBS



## Basic Information

Catalog Number:

85130-5-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11848

GenBank Accession Number:

BC050681

GeneID (NCBI):

10973

UNIPROT ID:

Q8N3C0

Full Name:

activating signal cointegrator 1  
complex subunit 3

Calculated MW:

111 aa, 13 kDa, 251 kDa

Observed MW:

245 kDa

Purification Method:

Protein A purification

CloneNo.:

242755D9

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

## Background Information

ASCC3 (Activating Signal Cointegrator 1 Complex Subunit 3) is a multifunctional protein involved in various cellular processes, including transcriptional regulation, DNA repair, and ribosome quality control. It is a member of the ASCC (Activating Signal Cointegrator 1 Complex) family, which functions as a bridge between co-repressors and co-activators in transcriptional regulation.

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85130-5-RR (ASCC3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85130-5-PBS in a different storage buffer formulation.