

For Research Use Only

# KCND2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85133-4-PBS



## Basic Information

Catalog Number:

85133-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15879

GenBank Accession Number:

BC110449

GeneID (NCBI):

3751

UNIPROT ID:

Q9NZV8

Full Name:

potassium voltage-gated channel,  
Shal-related subfamily, member 2

Calculated MW:

630 aa, 71 kDa

Purification Method:

Protein A purification

CloneNo.:

242754G12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01834-2, Kv4.2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85133-4-PBS. Detection antibody: 85133-3-PBS. Standard: Ag15879. Range: 1.563-100 ng/mL.