

For Research Use Only

TIMM44 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85134-1-PBS



Basic Information

Catalog Number:

85134-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4834

GenBank Accession Number:

BC033628

GeneID (NCBI):

10469

UNIPROT ID:

O43615

Full Name:

translocase of inner mitochondrial
membrane 44 homolog (yeast)

Calculated MW:

452 aa, 51 kDa

Purification Method:

Protein A purification

CloneNo.:

242774D12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01826-1, TIMM44 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85134-1-PBS. Detection antibody: 85134-3-PBS. Standard: Ag4834. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01826-2, TIMM44 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85134-1-PBS. Detection antibody: 85134-2-PBS. Standard: Ag4834. Range: 0.313-40 ng/mL