

For Research Use Only

# KIR2DL4/CD158d Recombinant antibody, PBS Only (Detector)

Catalog Number: 85135-3-PBS



## Basic Information

Catalog Number:

85135-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC041611

GeneID (NCBI):

3805

UNIPROT ID:

Q99706

Full Name:

killer cell immunoglobulin-like receptor, two domains, long cytoplasmic tail, 4

Calculated MW:

377 aa, 42 kDa

Purification Method:

Protein A purification

CloneNo.:

242235G1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01836-2, KIR2DL4/CD158d Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85135-1-PBS. Detection antibody: 85135-3-PBS. Standard: Eg1895. Range: 2.5-80 ng/mL.