

For Research Use Only

ACOX2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85136-2-PBS



Basic Information

Catalog Number:

85136-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11778

GenBank Accession Number:

BC047700

GeneID (NCBI):

8309

UNIPROT ID:

Q99424

Full Name:

acyl-Coenzyme A oxidase 2, branched chain

Calculated MW:

681 aa, 77 kDa

Purification Method:

Protein A purification

CloneNo.:

242752D11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01838-2, ACOX2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85136-2-PBS. Detection antibody: 85136-3-PBS. Standard: Ag11778. Range: 3.125-100 ng/mL.



Cytometric bead array standard curve of MP01838-1, ACOX2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85136-2-PBS. Detection antibody: 85136-1-PBS. Standard: Ag11778. Range: 1.563-100 ng/mL.