

For Research Use Only

# RNF128 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85137-2-PBS



## Basic Information

Catalog Number:

85137-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22834

GenBank Accession Number:

BC063404

GeneID (NCBI):

79589

UNIPROT ID:

Q8TEB7

Full Name:

ring finger protein 128

Calculated MW:

428 aa, 47 kDa

Purification Method:

Protein A purification

CloneNo.:

242759G3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926


E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.



# Selected Validation Data



Cytometric bead array standard curve of MP01819-2, RNF128 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85137-3-PBS. Detection antibody: 85137-2-PBS. Standard: Ag22834. Range: 0.625-40 ng/mL