

For Research Use Only

# TANK Recombinant antibody, PBS Only (Capture)

Catalog Number: 85138-1-PBS



## Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
85138-1-PBS	BC003388	Protein A purification
Concentration:	GeneID (NCBI):	CloneNo.:
1 mg/ml	10010	242767B11
Source:	UNIPROT ID:	
Rabbit	Q92844	
Isotype:	Full Name:	
IgG	TRAF family member-associated	
Immunogen Catalog Number:	NFκB activator	
AG25740	Calculated MW:	
	14 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01820-2, TANK Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85138-1-PBS. Detection antibody: 85138-2-PBS. Standard: Ag25740. Range: 0.625-40 ng/mL



Cytometric bead array standard curve of MP01820-1, TANK Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85138-1-PBS. Detection antibody: 85138-3-PBS. Standard: Ag25740. Range: 0.625-40 ng/mL