

For Research Use Only

TANK Recombinant antibody, PBS Only (Detector)

Catalog Number: 85138-2-PBS



Basic Information

Catalog Number:

85138-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25740

GenBank Accession Number:

BC003388

GeneID (NCBI):

10010

UNIPROT ID:

Q92844

Full Name:

TRAF family member-associated

NFkB activator

Calculated MW:

14 kDa

Purification Method:

Protein A purification

CloneNo.:

242767B12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01820-2, TANK Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85138-1-PBS. Detection antibody: 85138-2-PBS. Standard: Ag25740. Range: 0.625-40 ng/mL