

For Research Use Only

EPHB4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85139-2-PBS



Basic Information

Catalog Number:

85139-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10033

GenBank Accession Number:

BC052804

GeneID (NCBI):

2050

UNIPROT ID:

P54760

Full Name:

EPH receptor B4

Calculated MW:

987 aa, 108 kDa

Purification Method:

Protein A purification

CloneNo.:

242762C4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01825-1, EPHB4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85139-2-PBS. Detection antibody: 85139-1-PBS. Standard: Ag10033. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01825-2, EPHB4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85139-2-PBS. Detection antibody: 85139-3-PBS. Standard: Ag10033. Range: 0.781-100 ng/mL.