

For Research Use Only

# C1orf77 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85147-2-PBS



## Basic Information

Catalog Number:

85147-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26580

GenBank Accession Number:

BC002733

GeneID (NCBI):

26097

UNIPROT ID:

Q9Y3Y2

Full Name:

chromosome 1 open reading frame 77

Calculated MW:

248 aa, 26 kDa

Purification Method:

Protein A purification

CloneNo.:

242770E3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01863-2, C1orf77 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85147-1-PBS. Detection antibody: 85147-2-PBS. Standard: Ag26580. Range: 1.563-200 ng/mL.