

For Research Use Only

C1orf77 Recombinant antibody

Catalog Number: 85147-3-RR



Basic Information

Catalog Number:

85147-3-RR

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26580

GenBank Accession Number:

BC002733

GeneID (NCBI):

26097

UNIPROT ID:

Q9Y3Y2

Full Name:

chromosome 1 open reading frame 77

Calculated MW:

248 aa, 26 kDa

Observed MW:

26-30 kDa

Purification Method:

Protein A purification

CloneNo.:

242770H3

Recommended Dilutions:

WB 1:2000-1:10000

IF/ICC 1:500-1:2000

Applications

Tested Applications:

WB, IF/ICC, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, Jurkat cells, HEK-293T cells, THP-1 cells

IF/ICC : SH-SY5Y cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85147-3-RR (C1orf77 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (4% PFA) fixed SH-SY5Y cells using C1orf77 antibody (85147-3-RR, Clone: 242770H3) at dilution of 1:1000 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red).