

For Research Use Only

CRTC2,TORC2 Recombinant antibody, PBS Only (Detector)

Catalog Number:85150-3-PBS



Basic Information

Catalog Number:

85150-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3167

GenBank Accession Number:

BC053562

GeneID (NCBI):

200186

UNIPROT ID:

Q53ETO

Full Name:

CREB regulated transcription
coactivator 2

Calculated MW:

693 aa, 73 kDa

Purification Method:

Protein A purification

CloneNo.:

242808H1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01870-2, TORC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85150-1-PBS. Detection antibody: 85150-3-PBS. Standard: Ag3167. Range: 0.625-80 ng/mL