

For Research Use Only

KCNV1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85153-1-PBS



Basic Information

Catalog Number:

85153-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3376

GenBank Accession Number:

BC028739

GeneID (NCBI):

27012

UNIPROT ID:

Q6PIU1

Full Name:

potassium channel, subfamily V,
member 1

Calculated MW:

500 aa, 56 kDa

Purification Method:

Protein A purification

CloneNo.:

242865A2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01872-2, KCNV1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85153-3-PBS. Detection antibody: 85153-1-PBS. Standard: Ag3376. Range: 1.563-200 ng/mL.