

For Research Use Only

# ZFP36L2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85154-3-PBS



## Basic Information

Catalog Number:

85154-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13515

GenBank Accession Number:

BC005010

GeneID (NCBI):

678

UNIPROT ID:

P47974

Full Name:

zinc finger protein 36, C3H type-like 2

Calculated MW:

51 kDa

Purification Method:

Protein A purification

CloneNo.:

242776E7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01862-1, ZFP36L2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85154-3-PBS. Detection antibody: 85154-1-PBS. Standard: Ag13515. Range: 0.781-100 ng/mL