

For Research Use Only

NPM3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85156-3-PBS



Basic Information

Catalog Number:

85156-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2564

GenBank Accession Number:

BC054868

GeneID (NCBI):

10360

UNIPROT ID:

O75607

Full Name:

nucleophosmin/nucleoplasmin, 3

Calculated MW:

178 aa, 19 kDa

Purification Method:

Protein A purification

CloneNo.:

242863G2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01880-2, NPM3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85156-1-PBS. Detection antibody: 85156-3-PBS. Standard: Ag2564. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01880-1, NPM3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85156-2-PBS. Detection antibody: 85156-3-PBS. Standard: Ag2564. Range: 0.625-80 ng/mL