

For Research Use Only

SORBS1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85158-1-PBS



Basic Information

Catalog Number:

85158-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4828

GenBank Accession Number:

BC042612

GeneID (NCBI):

10580

UNIPROT ID:

Q9BX66

Full Name:

sorbin and SH3 domain containing 1

Calculated MW:

1292 aa, 142 kDa

Purification Method:

Protein A purification

CloneNo.:

242797A1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01871-2, SORBS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85158-2-PBS. Detection antibody: 85158-1-PBS. Standard: Ag4828. Range: 1.563-100 ng/mL.