

For Research Use Only

Mouse Cd44 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 85163-8-PBS



Basic Information

Catalog Number:

85163-8-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3159

GenBank Accession Number:

NM_009851.2

GeneID (NCBI):

12505

UNIPROT ID:

NP_033981

Full Name:

CD44 antigen

Calculated MW:

86 kDa

Purification Method:

Protein A purification

CloneNo.:

242666A8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CD44 concentration was determined to be 89.5 ng/mL in mouse serum.



Cytometric bead array standard curve of MP01873-4, MOUSE CD44 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85163-8-PBS. Detection antibody: 85163-9-PBS. Standard: Eg3159. Range: 0.781-100 ng/mL.