

For Research Use Only

# CD336 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85164-2-PBS



## Basic Information

Catalog Number:

85164-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_004828

GeneID (NCBI):

9436

UNIPROT ID:

O95944

Full Name:

natural cytotoxicity triggering  
receptor 2

Calculated MW:

31 kDa

Purification Method:

Protein A purification

CloneNo.:

242228C8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01875-2, CD336 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85164-2-PBS. Detection antibody: 85164-3-PBS. Standard: Eg1594. Range: 0.313-40 ng/mL