

For Research Use Only

SCARB1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85165-2-PBS



Basic Information

| | | |
|--------------------------------------|--|--|
| Catalog Number: 85165-2-PBS | GenBank Accession Number: BC022087 | Purification Method: Protein A purification |
| Concentration: 1 mg/ml | GeneID (NCBI): 949 | CloneNo.: 242842D4 |
| Source: Rabbit | UNIPROT ID: Q8WTV0 | |
| Isotype: IgG | Full Name: scavenger receptor class B, member 1 | |
| Immunogen Catalog Number: AG15759 | Calculated MW: 552 aa, 61 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01869-1, SR-BI Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85165-2-PBS. Detection antibody: 85165-1-PBS. Standard: Ag15759. Range: 1.563-200 ng/mL.