

For Research Use Only

TRIM8 Recombinant antibody, PBS Only

Catalog Number: 85169-1-PBS



Basic Information

Catalog Number:

85169-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26517

GenBank Accession Number:

BC021925

GeneID (NCBI):

81603

UNIPROT ID:

Q9BZR9

Full Name:

tripartite motif-containing 8

Calculated MW:

61 kDa

Observed MW:

61 kDa

Purification Method:

Protein A purification

CloneNo.:

242501B6

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

TRIM8 (Tripartite motif containing protein 8), also known as RNF27 (RING finger protein 27), a member of TRIM family proteins, is known to play a dual role as both tumor suppressor and oncogene, and to function at the crosstalk of cancer and innate immunity. TRIM8 can function at two subcellular sites-nucleus and cytoplasm-to regulate NF- κ B, one of the central signaling pathways that plays a critical role in carcinogenesis and inflammatory diseases (PMID: 33230447, 22262183). TRIM8 is composed of 551 amino acids (aa) with a molecular weight of 61.5 kDa, contains a Nuclear Localization Signal (NLS), and forms specific nuclear structures similar to the TRIM19/PML Nuclear Bodies (NBs) (PMID: 33807506, 29182544).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85169-1-RR (TRIM8 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85169-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 85169-1-RR against Human TRIM8 were performed. The affinity constant is 0.647 nM.