

For Research Use Only

GGGS Linker Recombinant antibody, PBS Only

Catalog Number: 85170-1-PBS



Basic Information

Catalog Number:

85170-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

242306C8

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

recombinant

Background Information

GGGS is one of the best fold-breaking linkers to expose fusion protein partners. This antibody detects the protein that contains GGGS linker.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human SCFV Lv-Hv Transfected HEK-293T cells were subjected to SDS PAGE followed by western blot with 85170-1-RR (GGGGS Linker antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85170-1-PBS in a different storage buffer formulation.

Immunohistochemical analysis of paraffin-embedded Transfected HEK-293T cells slide using 85170-1-RR (GGGGS Linker antibody) at dilution of 1:80000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 85170-1-PBS in a different storage buffer formulation.