

For Research Use Only

# Mouse CLEC12A Recombinant antibody, PBS Only (Capture)

Catalog Number: 85176-1-PBS



## Basic Information

Catalog Number:

85176-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2726

GenBank Accession Number:

/

GeneID (NCBI):

232413

UNIPROT ID:

Q504P2

Full Name:

C-type lectin domain family 12,  
member a

Calculated MW:

31 kDa

Purification Method:

Protein A purification

CloneNo.:

242620C12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01887-2, MOUSE CLEC12A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85176-1-PBS. Detection antibody: 85176-2-PBS. Standard: Eg2726. Range: 0.313-40 ng/mL