

For Research Use Only

Mouse Clec12a Recombinant antibody, PBS Only (Detector)

Catalog Number: 85176-2-PBS



Basic Information

Catalog Number:

85176-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

/

GeneID (NCBI):

232413

UNIPROT ID:

Q504P2

Full Name:

C-type lectin domain family 12,
member a

Calculated MW:

31 kDa

Purification Method:

Protein A purification

CloneNo.:

242620D1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01887-1, MOUSE CLEC12A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85176-3-PBS. Detection antibody: 85176-2-PBS. Standard: Eg2726. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01887-2, MOUSE CLEC12A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85176-1-PBS. Detection antibody: 85176-2-PBS. Standard: Eg2726. Range: 0.313-40 ng/mL