

For Research Use Only

HUWE1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85177-2-PBS



Basic Information

Catalog Number:

85177-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13763

GenBank Accession Number:

BC002602

GeneID (NCBI):

10075

UNIPROT ID:

Q7Z6Z7

Full Name:

HECT, UBA and WWE domain containing 1

Calculated MW:

437 aa, 482 kDa

Purification Method:

Protein A purification

CloneNo.:

242702E5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01895-1, HUWE1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85177-1-PBS. Detection antibody: 85177-2-PBS. Standard: Ag13763. Range: 0.625-80 ng/mL