

For Research Use Only

HUWE1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85177-4-PBS



Basic Information

Catalog Number:

85177-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13763

GenBank Accession Number:

BC002602

GeneID (NCBI):

10075

UNIPROT ID:

Q7Z6Z7

Full Name:

HECT, UBA and WWE domain
containing 1

Calculated MW:

437 aa, 482 kDa

Purification Method:

Protein A purification

CloneNo.:

242702A10

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01895-3, Human HUWE1 Recombinant Matched Antibody Pair - PBS only. 85177-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag13763. 85177-4-PBS was HRP conjugated as the detection antibody. Range: 1.25-80 ng/mL



The mean HUWE1 concentration was determined to be 36.46 ng/mL in HepG2 cell extract based on a 1.30 mg/mL extract load.