

For Research Use Only

# HOXA9 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85178-5-PBS



## Basic Information

Catalog Number:

85178-5-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13350

GenBank Accession Number:

BC010023

GeneID (NCBI):

3205

UNIPROT ID:

P31269

Full Name:

homeobox A9

Calculated MW:

30 kDa

Purification Method:

Protein A purification

CloneNo.:

242733B4

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP01901-3, Human HOXA9 Recombinant Matched Antibody Pair - PBS only. 85178-6-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag13350. 85178-5-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL.



The mean HOXA9 concentration was determined to be 31.36 ng/mL in HEK-293 cell extract based on a 4.50 mg/mL extract load.